


# RapidAI security and compliance

Providing AI-enhanced imaging with the most secure data management standards



RapidAI is committed to upholding the most rigorous information security standards and continually invests in people, processes and technology to protect the platform and your data. The Rapid platform has been developed using the National Institute of Standards and Technology (NIST) cybersecurity framework. It is also HIPAA compliant, GDPR compliant, ISO 13485, ISO 27001, ISO 27701 and SOC2 Type 1 certified.

The following safeguards are in place across the entire portfolio of products, ensuring international data and compliance requirements are being met at every level.

## Data security and compliance

- Data encrypted at rest using AES-256 encryption and in transit using TLS
- Immutable Operating System to strengthen the system against malware
- Containerized platform architecture to minimize attack surface
- Real-Time monitoring and logging with built-in telemetry
- Faster and simpler patching of application and security vulnerabilities
- Comprehensive Incident Management Plan and communication protocols
- Secure development practices including frequent vulnerability assessment of SBOM (software Bill of Materials)
- Role-based access control

## Mobile security

- Multi-factor authentication
- Biometric support
- Inactivity timeouts
- Integration with Mobile Device Management systems

## Information and utilization management

- User access and management controls

